





Large Area Infrared Systems

Improve quality of life for people with hearing loss using Contacta IR, providing clear sound, confidentiality and accessibility



Contents

	Section	Page
D _{IR}	Large Area Infrared Systems	
	Introduction	3
	Transmitters	4
	Receivers	8
	Charging	9
	Headphones	10
	Inductive Neck Loop	11



Infrared assistive listening systems provide clear, wireless hearing in large areas

Assistive technology systems provide independence, accessibility and well-being to users who require them.

Contacta IR assistive listening systems are ideal for when a sound source needs to be transmitted to multiple receivers, when ambient noise needs to be reduced, or when confidentiality is required within a confined space.

An audio source is converted to infrared light and projected to users via an **infrared transmitter**. **Infrared receivers** then convert the light back into audio, providing them with a simple and crystal clear listening experience.



Medium Area Infrared Transmitter with Modulator

The Medium Area Infrared Transmitter with Modulator combines an infrared modulator and transmitter, providing a wide angle of emission with a sleek and compact design (available in both black and white). With a coverage area of 180m², it is ideal for use in medium-sized environments.

All audio signals are modulated to infrared signals 2.3MHz or 2.8MHz.

The unit's range can also be increased and the angle enlarged by combining two Medium Area Infrared Transmitters. The range of the combined units is up to 370m² at one-channel transmission, ideal for large areas.

A wall and ceiling mounting kit is supplied.

Applications:

For use in any situation where infrared systems can provide multiple listeners with clear sound:

- Conference and Meeting Rooms
- · Cinemas and Theatres
- Courts
- Lecture Halls, Schools and Universities
- · Places of Worship

Colour Options

 The fascia is removable and available in both black and white

What's Included:

- Medium Area Infrared Transmitter w/ Modulator
- · Wall/Ceiling Mounting Kit

Accessories

• Portable Infrared Receiver (see page 8)

Technical Data:

Power Supply Input	100-240V AC / 50-60Hz
Power Supply Output	24V DC / 0.8A (2.5mm barrel plug - centre contact »positive«)
Connection Cable	2 x barrel plug, 2.5mm
Carrier Frequencies	2.3 MHz (A-channel) and 2.8 MHz (B-channel)
Transmission Power	0.7W
G	Single transmitter: 1-channel operation 180m²
Coverage Area	Combined mode (2 transmitters): 370m²
Line-Input	Cinch (RCA) for left and right and/or A and B-channel
Microphone Input	3.5 mm jack plug - sensitivity adjustable via control dial
Inputs/Outputs for	Primary Transmitter: Power-out 24V, audio-out "sync out"
Combined Mode	Secondary Transmitter: Power-in 24V, audio-in "sync-in"
LED Display	Input-level indicator LEDs, 1 per channel
Ambient Conditions	0 – +50°C ambient temperature (non-condensing, non-corrosive environment)
Mounting Options	Wall or ceiling mounting kit included

Physical Data:

•	
Dimensions	Height – 64mm (2.51") Width – 254mm (10") Depth – 79mm (3.11")
Weight	300g (0.66lbs)
Entre	White
Finish Bla	Black

Medium Area IR Transmitter with Modulator	IR-TX2
--	--------



Large Area Infrared Transmitter with Modulator

Our two-channel Large Area Infrared Transmitter with Modulator combines an infrared modulator and transmitter, generating a wide exit angle for the infrared signal while providing a sleek and compact design that eliminates the need for rack space. It is ideal for use in large environments.

Transmissions are at a bandwidth of 2.3MHz and 2.8MHz. In 1-channel mode, the transmitter supplies areas of up to $900m^2$. Additional transmitters allow the receiving range to be increased.

A wall and ceiling mounting kit is included.

Applications:

For use in any situation where infrared systems can provide multiple listeners with clear sound:

- Conference and Meeting Rooms
- Cinemas and Theatres
- Courts
- Lecture Halls, Schools and Universities
- Places of Worship

What's Included:

- Large Area Infrared Transmitter w/ Modulator
- Wall/Ceiling Mounting Kit

Accessories

• Portable Infrared Receiver (see page 8)

Technical Data:

Power Supply	Primary: 210 - 240V~ 50-60Hz Secondary: 28V DC
DC Input	Low-voltage connector, 28V DC, centre contact »negative«
Power Consumption	700mA
Carrier Frequencies	2.3 MHz (A-channel) and 2.8 MHz (B-channel)
Coverage Area	900m² (closed room)
Audio Inputs	Cinch (RCA) for left and right and/or A and B-channel
Microphone Input	6.3mm stereo-jack plug
	Audio-In: 18kΩ
Input Impedance	Mic-In: 1.3kΩ
Trigger Voltage of Audio-In	110mV
Input Level Display Mic-In	1.5mV
Inputs/Outputs for Combined Mode	Sync-in/Sync-Out Cinch (RCA) for left & right and/or A and B-channel
Display	Input level-LEDs, one per channel.
Ambient Conditions	0 – +40 °C ambient temperature (non-condensing, non-corrosive environment)
Mounting Types	Wall or ceiling mounting kit included

Physical Data:

Dimensions	Height – 327mm (12.8") Width – 412mm (16.2") Depth – 75mm (2.95")
Weight	1.1kg (2.42lbs)
Finish	Black

Large Area IR Transmitter with Modulator	IR-TX4
---	--------



Portable Infrared Receiver

Our Portable Infrared Receiver is a four channel receiver, compatible with our infrared transmitters using frequencies 2.3MHz and 2.8MHz. It features automatic gain control and high sensitivity for incoming infrared signals.

The receiver uses rechargeable batteries, and is a versatile and sturdy unit which requires either headphones or an assistive hearing neck loop.

Volume is amplified, adjustable up to 120dB to assist users with hearing loss.

What's Included:

- Portable Infrared Receiver
- Lanyard with Clip

Physical Data:

- Belt Clip
- 2x Rechargeable Batteries

Accessories

- 5 Bay Infrared Receiver Charger (see page 9)
- Inductive Neck Loop (see page 11)
- Headphones (see page 10)

Technical Data:

Receiver Frequencies	2.3MHz / 2.8MHz / 3.3MHz / 3.8MHz
Modulation	FM broadband
Battery Type	2x AAA rechargeable batteries; type NiMH, 1.2V 750 mAh
Frequency Response	60Hz to 12kHz ± 3 dB
Signal-to-Noise Ratio	54dB not weighted
Usage Period	20 hours
Maximum Output	6mW
C-#i	Volume Adjustment
Settings	4 channel selection - rotary control

Dimensions	Height – 112mm (4.40") Width – 70mm (2.75") Depth – 22mm (0.86")
Weight	200g (0.44lbs)
Finish	Black

Portable IR Receiver	IR-RX2
----------------------	--------

Single Bay Infrared Receiver Charger

IR-RX2-DC5

Our single bay receiver charger is available for our Portable Infrared Receiver (IR-RX2, see page 8).



What's Included:

- · Single Bay Infrared Receiver Charger
- Power Supply

Physical Data:

Dimensions	Height – 34mm (1.34") Width – 136mm (5.35") Depth – 85mm (3.35")
Weight	94g (0.21lbs)
Finish	Black

Technical Data:

Power Supply	Primary 100 -240 V DC 50-60 Hz / Secondary 12 V, 0.5 A, 6.0 W
--------------	---

Product Code:

Single Bay Infrared Receiver Charger	IR-RX2-DC1
--------------------------------------	------------

5 Bay Infrared Receiver Charger

IR-RX2-DC5

Our 5 Bay Infrared Receiver Chargers are available for the Portable Infrared Receiver (IR-RX2, see page 8).

This unit can hold up to 5 receivers in a standing position. Two chargers can be daisy-chained together to charge 10 receivers using a single power supply.



What's Included:

- 5 Bay Infrared Receiver Charger
- Power Supply

Physical Data:

Dimensions	Height – 38mm (1.49") Width – 160mm (6.29") Depth – 305mm (12.0")
Weight	540g (1.19lbs)
Finish	Black

Technical Data:

1 ever supply	Power Supply 12\	/ DC 1000mA
---------------	------------------	-------------

5 Bay Portable IR Receiver Charger	IR-RX2-DC5



Headphones

IR-HP1

Our black stereo headphones use a 2 metre lead and uses a standard 3.5mm audio jack.

Physical Data:

Dimensions	Height – 175mm (6.80") Width – 150mm (5.90") Depth – 15mm (0.59") Lead Length – 2 metres (78.7")
Weight	48g (0.10lbs)
Finish	Black

Technical Data:

Impedance	32Ω
Frequency Response	20-20000Hz
Sensitivity	104dB
Connection	3.5mm audio plug

Headphones	IR-HP1



Inductive Neck Loop

IR-NL1

The Inductive Neck Loop is a compact listening device compatible with our Portable Infrared Receiver (IR-RX2, see page 8). When connected to the receiver, the IR-NL1 will produce a magnetic field that lets users with Telecoil-enabled hearing instruments receive the sound directly in their hearing aids.

Physical Data:

Dimensions	Loop Cable Length – 800mm (31.5") Loop Cable Diameter – 4mm (0.15") Cord Cable Length – 450mm (17.7") Cord Cable Diameter – 3mm (0.12")
Weight	57g (0.12lbs)
Finish	Black

Technical Data:

Connection	3.5mm plug (stereo)
Input Power	250mW

Infrared Neck Loop	IR-NL1



Contact us now on:

\(\) +44 (0) 1732 223900 ⊠ sales@contacta.co.uk **\(\)** +44 (0) 1732 223909

Contacta Systems Ltd Office 13, Dana Estate Transfesa Road Paddock Wood, Kent TN12 6UT United Kingdom

www.contacta.co.uk